

Steritec Perfumed Virucidal Cleaner



Effective against
bacteria, yeasts
and enveloped
viruses in 60
seconds.

Effective against
non-enveloped
viruses in 60
minutes.

Fresh eucalyptus
fragrance, helps to
remove malodours,
leaving a clean,
fresh scent.



2 x 5L | STE-PVC-CON-2X5

A highly concentrated, perfumed disinfectant, suitable for use on most non-porous hard surfaces, leaving a clean, fresh scent after use. Giving excellent broad-spectrum activity against bacteria, yeasts, fungi and viruses. Ideal for use in hospitals, care homes where high levels of disinfection are required, as well as areas such as schools and leisure centres.

Not suitable for use in food preparation areas.



6 x 750ml HDR |
STE-PVC-RF-6X750

Directions for use

For professional use only.

1. Ensure that the surface is clean and free of debris before application.
2. Always dilute before use.
3. For efficacy against Enveloped Viruses and Norovirus, dilute 30ml of product with clean water into a 750ml trigger spray (4%).
Additional dilution rates and contact times can be found on the next page, depending on efficacy required.
4. Apply to surface to be cleaned via spray or by wiping, paying particular attention to corners and under edges.
5. Ensure that the surface is completely wet after application.
6. Leave for up to 5 minutes (4% dilution) before rinsing or allow to air dry.
Surfaces that come into contact with food must be rinsed after application.

Product Info

Colour: **Turquoise**

Fragrance: **Eucalyptus**

pH: **12.5 – 13.5**

Shelf Life: **36 Months**

CLEENOL

Tel: +44(0) 1295 251 721

Email: sales@cleenol.co.uk

Neville House, Beaumont Road, Banbury, Oxon. OX16 1RB.

Cleenol Group Ltd is a limited company registered in England and Wales under CRN 00635803.

Efficacy Summary

Bacteria

EN TEST	DILUTION (%)	CONTACT TIME	CONDITION	ORGANISM	
EN 1276	2%	60 seconds	Dirty	Pseudomonas aeruginosa Staphylococcus aureus	Escherichia coli Enterococcus hirae
EN 1276	0.5%	5 mins	Dirty	Listeria monocytogenes	Salmonella typhimurium
EN 13727	2%	60 seconds	Dirty	Pseudomonas aeruginosa Staphylococcus aureus	Escherichia coli Enterococcus hirae
EN 13697	2%	5 mins	Dirty	Pseudomonas aeruginosa Staphylococcus aureus	Escherichia coli Enterococcus hirae
EN 16615	1%	60 seconds	Dirty	Pseudomonas aeruginosa Staphylococcus aureus	Enterococcus hirae Proteus hauseri
EN 14561	2.5%	15 mins	Clean	Pseudomonas aeruginosa Staphylococcus aureus	Enterococcus hirae

Viruses

EN TEST	DILUTION (%)	CONTACT TIME	CONDITION	ORGANISM	
EN 14476	4%	60 seconds	Clean	Modified Vaccinia virus ankara	Bovine Corona virus
EN 14476	3%	5 mins	Clean	Murine Norovirus Influenza A (H1N1/H7N9)	Adenovirus
EN 14476	5%	5 mins	Dirty	Murine Norovirus	Adenovirus
EN 14476	5%	60 mins	Dirty	Poliovirus	
EN 16777	4%	60 seconds	Clean	Modified Vaccinia virus ankara	
EN 16777	5%	30 mins	Clean	Murine Norovirus	Adenovirus
EN 17111	1%	15 mins	Clean	Modified Vaccinia virus ankara	

Fungi

EN TEST	DILUTION (%)	CONTACT TIME	CONDITION	ORGANISM	
EN 16615	2%	15 mins	Clean	Aspergillus brasiliensis	
EN 1650	4%	30 mins	Clean	Aspergillus brasiliensis	

Yeasts

EN TEST	DILUTION (%)	CONTACT TIME	CONDITION	ORGANISM	
EN 1650	2%	60 seconds	Dirty	Candida albicans	
EN 13624	4%	60 seconds	Dirty	Candida albicans	
EN 13697	3%	5 mins	Dirty	Candida albicans	
EN 16615	1%	60 seconds	Dirty	Candida albicans	
EN14562	3%	15 mins	Clean	Candida albicans	

Use biocides safely and sustainably. Always read the label and product information before use.

CLEENOL

Tel: +44(0) 1295 251 721

Email: sales@cleenol.co.uk

Neville House, Beaumont Road, Banbury, Oxon. OX16 1RB.

Cleenol Group Ltd is a limited company registered in England and Wales under CRN 00635803.