

Evolution EV3

Antibacterial Multi-Surface Cleaner

EVOLUTION
SUPER CONCENTRATES

Neutral pH, making it suitable for use on most hard surfaces.

Effective against bacteria in 30 seconds.

Fragrance-free formulation, ideal for use in food preparation areas.



2 x 1.5L | EV3

A highly concentrated and economical disinfectant cleaner. It's highly effective formulation allows for easy cleaning and disinfection of most non-porous hard surfaces. It's fragrance-free formulation makes it ideal for use in food preparation areas, including on worktops, microwaves, chopping boards and stainless steel surfaces.



6 x 750ml |

EVOESD-3-RF-6X750

Directions for use

For professional use only.

1. Always dilute before use.
2. Dispense 20ml product into a 750ml refill flask (2.7%) and fill with clean water.
3. Apply product directly to the surface or onto a cloth and wipe, paying particular attention to edges and corners.
4. Ensure that the surface is completely wet after application.
5. Leave for up to 30 seconds before rinsing or allow to air dry.
Surfaces that may come into contact with food must be rinsed after application.

Reloading Pouch Pack

1. Remove empty pouch from dispenser and dispose of responsibly.
2. Remove screw cap from new pouch.
3. Insert pouch into Evolution Dispenser, and depress locking cap firmly into pouch.
4. Lock cabinet. Pump until product is dispensed.

Product Info

Colour: **Red**

Fragrance: **None**

pH: **7 - 8**

Shelf Life: **24 Months**

CLEENOL

Tel: +44(0) 1295 251 721

Email: sales@cleenol.co.uk

Neville House, Beaumont Road, Banbury, Oxon. OX16 1RB.

Cleenol Group Ltd is a limited company registered in England and Wales under CRN 00635803.

Efficacy Summary

Bacteria

EN TEST	DILUTION (%)	CONTACT TIME	CONDITION	ORGANISM	
EN 1276	2.7%	30 seconds	Clean	Pseudomonas aeruginosa Staphylococcus aureus	Escherichia coli Enterococcus hirae

Use biocides safely and sustainably. Always read the label and product information before use