

edencleen All Purpose Sanitiser

Effective against
bacteria and
enveloped viruses
in 30 seconds.

Ready to use
product for ease of
use.

Fragrance-free
formulation, ideal
for use in food
preparation areas.

An effective ready-to-use disinfectant cleaner, for use on most non-porous hard surfaces, including varnished wooden tables. Ideal for use in areas where easy and efficient cleaning and sanitising is required. It's fragrance-free formulation makes it suitable for use in food preparation areas. Formulated with plant-based materials to increase sustainability of the product.

Supplied in a 5L refill pack to reduce plastic waste.

Directions for use

For professional use only.

1. Dispense product into the refill flask – do not dilute.
2. Always test an inconspicuous area before use.
3. Apply to surface to be cleaned via spray or by wiping.
Paying particular attention to corners and under edges.
4. Ensure the surface is completely wet after application for at least 30 seconds.
5. Rinse or allow to air dry.

eden
cleen



2 x 5L | EDN-APS-2X5



6 x 750ml HDR |
EDN-APS-RF-6X750

Product Info

Colour: **Colourless**

Fragrance: **None**

pH: **2 – 3**

Shelf Life: **24 Months**

CLEENOL

Tel: +44(0) 1295 251 721

Email: sales@cleenol.co.uk

Neville House, Beaumont Road, Banbury, Oxon. OX16 1RB.

Cleenol Group Ltd is a limited company registered in England and Wales under CRN 00635803.

Efficacy Summary

Bacteria

EN TEST	DILUTION (%)	CONTACT TIME	CONDITION	ORGANISM	
EN 1276	N/A	30 seconds	Clean	Pseudomonas aeruginosa Escherichia coli	Staphylococcus aureus Enterococcus hirae
EN 13697	N/A	5 mins	Clean	Pseudomonas aeruginosa Escherichia coli	Staphylococcus aureus Enterococcus hirae

Viruses

EN TEST	DILUTION (%)	CONTACT TIME	CONDITION	ORGANISM	
EN 14476	N/A	30 seconds	Clean	Vaccinia virus	

Use biocides safely and sustainably. Always read the label and product information before use.