

PRODUCT USE GUIDE FOR NOVACROSS BN115

Product Description

A multipurpose disinfectant cleanser and deodoriser. Cleans, disinfects and deodorises, all in one operation. It is based on quaternary ammonium compounds, cleaning chemicals and reodorisers. Designed for routine use in hospitals, nursing homes, clinics, surgeries, kitchens, canteens, hotels, schools and offices; for cleaning and disinfecting floors, walls, paintwork, equipment, sinks, drains, washrooms, lavatories, waste bins, etc.

Product Features

Bactericide

- Contains 12,500 ppm of a powerful quaternary ammonium compound germicide. Chosen for its high activity, the broad-spectrum bactericidal system has wide-range killing powers, effective against both Gram positive and Gram negative organisms. Penetrates cracks, crevices and porous surfaces.
- Novacross complies with the European Suspension Test BSEN1276 and is effective against methicillin resistant staphylococcus aureus (MRSA), salmonella typhimurium, listeria monocytogenes, pseudomonas aeruginosa, Escherichia coli, staphylococcus aureus and enterococcus hirae with a contact time of 5 minutes.

Powerful Cleanser

- Contains carefully balanced cleaning components avoiding the use of highly alkaline compounds, this neutral cleanser is quickly effective against a wide range of soils. Safe to use on all hard surfaces which water alone will not harm. Efficient in hot or cold, hard or soft water.

Deodoriser/Perfume

- Effective cleansing will inhibit or destroy odour-producing bacteria. In addition, an unobtrusive but refreshing reodoriser lingers to maintain a pleasant working atmosphere, the signal of a clean environment.
- The product is also available unperfumed, and also with a pleasant citrus lemon perfume.

Direction for Use

Dilution Rates

- | | | |
|--------------------------------|------------|-----------------------|
| • For Heavy Duty Cleaning | Use 1 : 20 | (500ml per 10 litres) |
| • For General Routine Cleaning | Use 1 : 40 | (250ml per 10 litres) |
| • For Light Duty Cleaning | Use 1 : 50 | (200ml per 10 litres) |



Cleenol Group Ltd

Beaumont Road, Banbury, Oxon, OX16 1RB. UK.
Tel: +44 (0) 1295 251721 Fax: +44 (0) 1295 269561
Email: sales@cleenol.co.uk
www.cleenol.co.uk

JH/SMW/ES



Certificate No. 316



Certificate No. E3539

PRODUCT USE GUIDE FOR NOVACROSS BN115

Sickrooms, Wards, Health Centres, Clinics or Nurseries

- For the washing and disinfecting of floors, walls, ceilings, paintwork, furniture, utensils and equipment. Also for damp dusting of furniture, ledges, rails, toys and equipment (wipe with a cloth wrung out in a solution of 'NOVACROSS').

Kitchens, Canteens, Hotels, etc.

- For the washing and disinfecting floors, walls, ceilings, paintwork, equipment, etc.

Toilets, Washrooms, etc.

- For the washing and disinfecting of floors, walls, doors, ceilings, paintwork, etc. Use for cleaning toilet seats, cisterns, fittings, handles, washbasins, showers, etc.
- ('NOVABLOO' Toilet Bowl Cleanser, an acidic disinfectant product, is specially recommended for the cleaning, descaling and disinfection of the actual WC bowl and for urinals).

For washbasins, etc.

- Use of a hand-spray and special scrubbing sponge will speed the operation considerably.

Flushing Sinks, Drains, Overflows, Waste Bins, etc.

- A more concentrated solution may be used for disinfecting the problem areas. May also be used in lavatories, urinals, etc.

For Use in Scrubber Driers

- Use rate: dilute 2% in warm water. Add to the detergent tank.

Disinfection of Mops, Cloths, Cleaning Equipment, etc.

- Soak and wash out mops, cleaning cloths, pads, brushes, etc. in a solution of 'NOVACROSS' to keep them hygienic and germ-free. Wipe mechanical equipment such as vacuum cleaners, cleaning machinery, buckets, etc. with a cloth wrung out in a 1-20 solution.

Handling

- As supplied, 'NOVACROSS' Sanitizing Fluid may be regarded as having a low order of irritation. Diluted, it can be regarded as harmless when used in a sensible and responsible manner.

Storage

- Store in cool conditions away from direct sources of heat
- Avoid extremes of temperature
- Do not allow to freeze
- Ensure the cap is tightly fitted after use
- Keep out of reach of children



Cleenol Group Ltd

Beaumont Road, Banbury, Oxon, OX16 1RB. UK.
Tel: +44 (0) 1295 251721 Fax: +44 (0) 1295 269561
Email: sales@cleenol.co.uk
www.cleenol.co.uk

JH/SMW/ES



Certificate No. 316



Certificate No. E3539

PRODUCT USE GUIDE FOR NOVACROSS BN115

Shelf Life

- Shelf life 24 months, if stored as recommended.
- Use within 6 months of opening.

Environmental

- Do not allow large quantities of neat liquid to enter surface waterways.
- This product uses recyclable packaging. The detergents used in the formulation are biodegradable to EU guidelines under the detergents regulations (EC) No 648/2004.
- It contains no Phosphates.
- None of the components of this product are animal derived.

For further information please scan your QR Smart Reader here



Code Nos.

BN1150	2x5L
BN1151/6	6x1L
BN115F	6x750ml refill flasks

Data Sheet No. BN115

Note: This product is only to be used for the purpose for which it was intended and in accordance with the specified instructions. This product should never be mixed with other chemical products and should be handled with care and stored carefully out of reach of children.



Cleenor Group Ltd

Beaumont Road, Banbury, Oxon, OX16 1RB. UK.
Tel: +44 (0) 1295 251721 Fax: +44 (0) 1295 269561
Email: sales@cleenol.co.uk
www.cleenol.co.uk

JH/SMW/ES



Certificate No. 316



Certificate No. E3539